

XpreciaPrime™ 4U

Coagulation analyzer for easy self-monitoring
anytime, anywhere

Technical specifications

Analyzer features

- 1 Test strip eject button
- 2 Touchscreen display
- 3 Power button
- 4 Barcode scanner
- 5 Test strip port
- 6 USB port



Affordable

- Avoid travel expenses and test whenever convenient
- Take control of your health without worrying about your budget



Accurate

- Correlation of results with laboratory analysis
- Reliable results, less repetition
- Automatic strip quality check



Easy, intuitive

- User friendly software and touchscreen
- Automatic calibration - no reagents needed



Rapid results

- Results obtained in under a minute
- No more waiting for lab results



Portable

- Fits easily in your hand
- Micro-USB connectivity
- Rechargeable battery



Safe to use

- Eject button for safe handling

Note: This product is not currently available for sale in the USA

Technical Specifications

Sample

Sample type	Droplet of blood, note, no lancing device is supplied with this product
Sample volume	8 µL
Results displayed	'International Normalized Ratio' (INR) and 'prothrombin time' (PT) in seconds

Haematocrit range 25 – 55% inclusive

Reporting range 0.8 – 8.0 INR

Features

Automatic power off	Configurable: disabled, 5 mins, 15 mins, 30 mins
Dimensions	147 × 84 × 32 mm (5.79 × 3.31 × 1.26 in)
Design lifetime	3 years

Weight 210 g (7.41 oz)

Memory recall 1000 test results

Analyzer/system operating conditions

Temperature	15°C – 32°C (59°F – 89°F)
Altitude	Electrical safety tested up to 2000 m

Relative humidity < 80% (without condensation)

Transport and storage conditions (in packaging)

Temperature	-20°C – 40°C (-4°F – 104°F)
Altitude	Electrical safety tested up to 2000 m

Relative humidity < 75% (without condensation)

Test strips storage conditions

Temperature	2°C – 30°C (35.6°F – 86°F)
--------------------	----------------------------

Relative humidity < 75% (without condensation)

Supplied accessories

1 x	Xprecia Prime™ 4U analyzer
1 x	Replacement cap
1 x	User guide

1 x	Charger
1 x	Power supply (including USB cable and 4 sockets)
1 x	Quick start guide

General

Input	5V === 0.9A
Capacity	Rechargeable 1.85Ah (LI-ION) with power monitoring
Pollution degree rating	Pollution degree II
Connectivity	Micro-USB

Rating 3.7V === 0.5A MAX

Water ingress rating IPX0

Overvoltage category Category I

Xprecia Prime™ and all related names, logos, product and service names, designs, and slogans are trademarks of Universal Biosensors Pty Ltd and may not be used without the prior written permission of Universal Biosensors Pty Ltd. All other names, logos, product and service names, designs, and slogans are the property of their respective owners.

Note: This product is not currently available for sale in the USA.

